

purify and capture antibody fragments with a new synthetic ligand affinity adsorbent

Fabsorbent™ F1P HF

In response to the growth of antibody fragments as research tools and biotherapeutic products, ProMetic BioSciences Ltd (PBL) has developed a new affinity adsorbent for the capture and purification of a wide range of antibody fragments.

Fabsorbent™ F1P HF is an adsorbent with a synthetic ligand which can be used as a superior alternative to Protein L for the capture and purification of antibody fragments including monovalent antibody fragments (e.g. Fab, scFv), engineered antibody variants and single domain antibodies.

The robustness and broad applicability of Fabsorbent™ F1P HF makes this adsorbent the ideal choice for both early stage research and development applications and also for larger-scale manufacturing uses where the full benefit of the high-flow/low pressure PuraBead® HF support matrix can be utilised. Given its broad selectivity, Fabsorbent™ F1P HF provides a “platform” technology for the purification of antibody fragments.

KEY BENEFITS

PROVEN TECHNOLOGY

Sharing the same chemistry as PBL's proven Mimetic Ligand™ adsorbents which are used extensively in commercial manufacturing processes, Fabsorbent™ F1P HF has been developed in close collaboration with antibody fragment producer companies. Using PBL's Chemical Combinatorial Library CCL® technology and extensive computer modelling and ligand screening expertise, PBL have identified the Fabsorbent™ F1P ligand which binds selectively to kappa (κ) and lambda (λ) light chains of a variety of mammalian antibodies including human, ovine, murine and bovine sources.

UNIVERSALLY APPLICABLE

Binds and purifies most antibody species (antibody fragments and whole antibodies) that contain κ and/or λ light chains.

ROBUST

Fabsorbent™ F1P HF is fully resistant to concentrated sodium hydroxide (up to 1 M) enabling effective column cleaning and sanitisation thereby prolonging column lifetime.

SAFE

Synthetic ligand not derived from cell culture. No components of animal origin are used in the manufacture of Fabsorbent™ F1P HF. This product is manufactured in PBL's ISO 9001 certified facility and is supported by a comprehensive Regulatory Support File.

OPERATIONAL FLEXIBILITY

PuraBead® HF base matrix enables the operation of large diameter columns at high linear flow rates (> 500 cm/h).

MULTIPLE FORMATS FOR EASE OF USE

Fabsorbent™ F1P HF is available in a variety of easy to use formats including 96-column block format PuraPlate™ adsorbent screening kits, pre-packed columns and a variety of slurry pack sizes.

ANTIBODY FRAGMENTS

- Fab' expressed in *E.coli*.
- Mammalian Fab fragments from papain catalysed antibody digests.
- Mammalian F(ab')₂ fragments from pepsin catalysed antibody digests.
- Human antibody λ and κ light chains.
- scFv expressed in *E.coli*.
- Single domain antibody fragments expressed in *E.coli*.
- Full length monoclonal antibodies.
- Human, bovine, ovine and murine antibodies and antibody fragments.

F(ab')₂

F(ab')₂ antibody fragments are prepared by pepsin cleavage of full-chain IgG which cleaves the Fc region to yield two Fab' fragments linked together. Fabsorbent™ F1P HF may be used to capture and purify F(ab')₂ fragments and, since pepsin is not bound by Fabsorbent™ F1P HF, this adsorbent may be used to simultaneously purify F(ab')₂ fragments and remove free pepsin (Figures 1 and 2).

FIGURE 1

Capture and elution of F(ab')₂ from pepsin digest.

Column volume: 10 mL; flow rate: 250 cm/h; equilibration buffer: 25 mM Tris, pH 8.0; elution buffer: 50 mM sodium citrate, pH 3.0; CIP: 0.5 M NaOH.





FIGURE 2

Non-reduced SDS-PAGE of F(ab')₂ from pepsin digest. Lane 1: molecular weight marker; Lane 2: IgG standard; Lane 3: pepsin digest; Lane 4: flow through and post load wash; Lane 5: elution.



DOMAIN ANTIBODY FROM *E. COLI* LYSATE

V_L domain antibodies are the smallest antigen binding antibody fragments. Recombinant V_L antibodies are commonly produced in microbes such as *E. coli* with the resulting fragment being purified from cell lysate. Fabsorbent™ F1P HF has been used to successfully capture and purify V_L domain antibodies (Figure 3). Non-reduced SDS-PAGE of the fractions from Fabsorbent™ F1P HF (Figure 4) demonstrates the capture and purification of V_L domain antibody from *E. coli* lysate.

FIGURE 3

Capture and elution of a domain antibody from *E. coli* lysate. Column volume: 1 mL; 3 minute residence time; equilibration buffer: 25 mM Tris, pH 9.0; elution buffer: 50 mM sodium citrate, pH 3.0; CIP: 0.5 M NaOH.





FIGURE 4

Non-reduced SDS-PAGE. Lane 1: molecular weight marker; Lane 2: V_L domain *E. coli* lysate; Lane 3: flow through; Lane 4: post load wash; Lane 5: elution; Lane 6: CIP.



MONOCLONAL IgG₁

Fabsorbent™ F1P HF was designed for the capture of antibody fragments and can therefore be used for the capture and purification of whole molecule IgG. Fabsorbent™ F1P HF can capture and purify IgG₁ directly from CHO cell culture supernatant (Figure 5 and 6). Fabsorbent™ F1P HF is designed for the capture of light chain containing fragments. All antibody related proteins in CHO cell culture supernatant will be captured and eluted by Fabsorbent™ F1P HF.

FIGURE 5

Chromatogram of the capture and elution of IgG₁ from CHO cell culture supernatant. Column volume: 2 mL; 3 minute residence time; equilibration buffer: PBS, pH 6.0; elution buffer: 50 mM sodium citrate, pH 3.0; CIP: 0.5 M NaOH.





FIGURE 6

Non-reduced SDS-PAGE; Lane 1: molecular weight marker; Lane 2: load; Lane 3: elution; Lane 4: CIP.



PROPERTIES

Ligand	F1P synthetic triazine ligand
Matrix	PuraBead® P6HF (6% near monodisperse highly cross-linked agarose).
Particle size	90 ± 10 µm
Binding capacity	Up to 60 mg IgG/mL and 20 mg Fab/mL resin (dependant on antibody/fragment type and source).
Operational Flow rate	Up to 500 cm/h
Operating pressure	Up to 3 bar (45 psi)
Operating pH	pH 1.5 - pH 14.0 (Intermittent)
pH stability	Long term (3 months) pH 3 - pH 13
Chemical stability	All commonly used aqueous buffers and co-solvents.
Sanitisation	0.5 - 1.0 M sodium hydroxide, 25°C
Storage	20% ethanol

ORDER

FABSORBENT™ F1P HF SLURRY FOR COLUMN PACKING

Fabsorbent™ F1P HF (25 mL)	3904-00025
Fabsorbent™ F1P HF (100 mL)	3904-00100
Fabsorbent™ F1P HF (500 mL)	3904-00500
Fabsorbent™ F1P HF (1000 mL)	3904-01000

PBL offers a range of larger pack sizes for supply of bulk resins into development and regulated cGMP manufacturing processes. For more information on Fabsorbent™ F1P HF and supply related matters please do not hesitate to contact us at: sales@prometicbiosciences.com

FABSORBENT™ F1P HF PRE-PACKED COLUMNS

(for attachment to chromatography workstations)

F1P HF Pre-Packed Columns (5 x 1 mL)	4904-00001
--------------------------------------	------------

FABSORBENT™ F1P HF PURAPLATE™

(96 column plate for rapid antibody screening and purification method development)

F1P HF PuraPlate™ (96 column plate)	2233
-------------------------------------	------

CONTACT US

At www.prometicbiosciences.com.

For sales questions, contact: sales@prometicbiosciences.com.

For tech support questions, contact: techsupport@prometicbiosciences.com.

TEL +44 (0) 1223 420300/FAX +44 (0) 1223 420270

Fabsorbent™, Mimetic Ligands™, PuraPlate™, ProMetic BioSciences Ltd and the PBL Logo are trademarks of ProMetic BioSciences Ltd.

Chemical Combinational Library CCL® is registered with the U.K. Patent & Trademark Office.

PuraBead® is registered with the U.S. Patent & Trademark Office.

© ProMetic BioSciences Ltd 2010. Issue - 060510